

Cat. No. B10-2638-W0 GENECON Due with Pulley INSTRUCTION MANUAL

Thank you for purchasing the GENECON (DC generator) DUE with a pulley. This GENECON has a pulley instead of a handle. It may be used for energy conversion experiments such as hoisting a weight hung from the pulley and, conversely, dropping the weight hung from the pulley to light a miniature lamp by connecting it to another GENECON.

[Instructions]

- · Do not bring the product close to fire or moisten it.
- · Do not store the product in a hot and humid place.
- · Be very careful of shocks when dropping a weight. If the shocks are too intense due to heavy weight, the GENECON may break.
- · Before starting an experiment, secure the GENECON DUE with a pulley to a stand firmly and make sure that the stand is also secured on flat surface.
- · Fix a weight firmly to the string.
- · Tighten the pulley firmly with a hexagon wrench if it becomes loose.
- If a weight is dropped from the pulley in the condition where no load is connected or a low-load LED is connected to the pulley, it drops very quickly. Be very careful.

Product Specifications

Materials: Body (Polycarbonate)

Pulley (Aluminum)

Sizes: $115 \times 140 \times 43 \text{ mm}$

 $\phi 50 \times 10 \text{ mm (pulley)}$

Accessory: Kite string, 2 m



[Usage]

- (i) Attach the kite string supplied with the product to the pulley. The pulley has a hole for knotting the string in the outer circumference. Pass the string through it and fasten it. Wind the string around the pulley before the experiment every time.
- (ii) Plug the attached cable to the body (red cap). Connect a miniature lamp, etc to it and you can start using it for experiments.



[Examples of Experiments]

Example 1: Experiment of hoisting a weight hung from the pulley

- 1. Connect a large weight (500g) or PET bottle filled with water to the top of the string attached to the pulley.
- 2. Connect the GENECON DUE with pulley with another GENECON DUE by its alligator clips using the attached cables.
- 3. Secure the GENECON DUE with pulley to a steel stand or hold firmly with a hand in the condition where the weight can be hoisted when the pulley rotates.
- 4. Rotate the handle of the connected GENECON DUE. Make sure that the pulley rotates, and the weight is wound up. It is possible to check the heaviness of the





Fig. 1

GENECON DUE handle or calculate the transition of energy if a voltmeter or ammeter is connected.

Example 2: Dropping the weight hung from the pulley to light a miniature lamp

- 1. Connect a miniature lamp to the cable connected to the GENECON DUE with pulley. Miniature bulb with specifications 3.8V, 0.3A or equivalent is recommended.
- 2. Attach a large weight of 1 kg or so or a PET bottle filled with water to the top of the string and support firmly.
- Secure the GENECON DUE with pulley with a steel stand, etc. or hold it with a hand firmly.
 When it is fixed with a steel stand, the stand may fall.
 Secure the base with tape, etc.
- 4. Wind the kite string with the weight around the pulley and drop the weight. The pulley rotates and the miniature lamp will light.
- 5. It is possible to calculate energy conversion if the dropping distance of 1 meter is measured, a voltmeter or ammeter is connected, and data are recorded.



* The dropping speed of the weight differs based on the specification of the miniature lamp.

NaRika Corporation

5-3-10 Sotokanda, Chiyoda, Tokyo, 101-0021 Japan http://www.global.narika.jp/